

Tesamorelin-Ipamorelin-Blend

PEPTIDES

DOSING CYCLING: 12-16 week cycle followed by 4-8 week washout period. Shorter cycles will not yield tangible results- do not begin this product with the intention of “trying it out” for 3 or 4 weeks. This will not be effective. Each week should consist of 5 consecutive dosing days

DOSING CYCLING: 12-16 week cycle followed by 4-8 week washout period. Shorter cycles will not yield tangible results- do not begin this product with the intention of “trying it out” for 3 or 4 weeks. This will not be effective. Each week should consist of 5 consecutive dosing days followed by 2 “rest” days.

TIMING: Dose at night, at least 90 minutes after your last meal. **RECONSTITUTION SOLUTION(S):** Hospira brand BAC water preferred (many companies do not maintain a known, consistent pH level like Hospira does) .6% Acetic Acid Solution Hospira brand BAC water is preferred because of the absolute consistency of its pH level- if the BAC water is a bit too acidic or too basic , peptides can lose their electrical charge, causing them to clump together. Extreme pH differences can denature the peptide, altering its structure and potentially rendering it ineffective. Acetic Acid solution is mandatory in order to prevent gelling of the product **ADMINISTRATION:** Subcutaneous Injection **STORAGE:** Store in refrigerator prior to reconstitution. When you are ready to prepare a new vial, remove the peptide from the refrigerator and allow it to come to room temperature before mixing with room temperature acetic acid solution and bacteriostatic water. After reconstitution, Tesamorelin/Ipamorelin blend is to be stored at room temperature: do not refrigerate . **POST-RECONSTITUTION EXPIRATION DATE:** 2 weeks.

RECONSTITUTION Mix with .5mL Acetic Acid until dissolved; then add 2.5mL Bacteriostatic Water

DOSING : TESAMORELIN 10MG + IPAMORELIN 3MG Weeks 1-6: Take 23 units once a day, at bedtime, 5 consecutive days a week Weeks 7-12: Either remain at 23 units per day/5 days per week for the duration of the cycle OR increase to 46 units per day/5 consecutive days each week for the duration of the cycle

RESEARCH USE DISCLAIMER: All information in this document is provided strictly for laboratory research and educational purposes. It is not intended for human or veterinary use, diagnosis, or treatment.